



the **SUB™ FLUSH KIT**  
*A Special Flush Kit with T-FloLoc™*  
*for use with the SUB System*



**A New Option for Preventing Occlusion and Biofilm Formation in the SUB System**  
*designed and developed in collaboration with veterinarians*

This convenient kit contains the items you need to quickly and effectively flush the SUB™ System. Conveniently packaged into one easy-to-use, compact kit, this set saves you precious time and packs the power of **T-FloLoc™**, the PREMIUM catheter flushing and locking solution offering both patency and infection control.

**SUB™ FLUSH KIT CONTENTS (5 PER CASE)**

**SFK-22** - for use with All Shunting Ports

- 1 x T-Port Connector
- 1 x 3-way Stop-Cock
- 1 x 22-Gauge Huber Point Needle
- 1 x 3mL Syringe
- 1 x 2.5mL Sterile Saline in 3mL Syringe
- 1 x 2mL **T-FloLoc™** in 12mL Syringe
- Instructions For Use (1 per Case of 5 Kits)

**SFK-20** - option for Swirl and Large Shunting Ports

- 1 x T-Port Connector
- 1 x 3-way Stop-Cock
- 1 x 20-Gauge Huber Point Needle
- 1 x 3mL Syringe
- 1 x 2.5mL Sterile Saline in 3mL Syringe
- 1 x 2mL **T-FloLoc™** in 12mL Syringe
- Instructions For Use (1 per Case of 5 Kits)



**ORDERING INFORMATION**

| Catalog Number        | Description                           | Price            |
|-----------------------|---------------------------------------|------------------|
| SFK-22-01 / SFK-20-01 | <u>1 Case</u> of 5 x SUB™ Flush Kits  | Call For Pricing |
| SFK-22-05 / SFK-20-05 | <u>5 Case</u> of 5 x SUB™ Flush Kits  | Call For Pricing |
| SFK-22-10 / SFK-20-10 | <u>10 Case</u> of 5 x SUB™ Flush Kits | Call For Pricing |

*NOTE: each flush kit is sterile packaged, suitable for sterile field use*



For Technical Information, please contact:

Allyson Berent, DVM, DACVIM  
 allyson.berent@gmail.com

Chick Weisse, VMD, DACVS  
 chick.weisse@gmail.com

**Norfolk Vet Products Inc.**

7350 North Ridgeway Avenue • Skokie, Illinois 60076 USA

Tel: 847.674.7143

Fax: 847.674.7066

E-mail: info@norfolkvetproducts.com